

SW6

VER:1.0

USER MANUAL



Warning

Dear user:

Thanks for using our company's products, in order to fully make use of the products, we sincerely suggests follow:

- ◆◆ Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆ Please follow the instruction to fit and operate the machine correctly.
- ◆◆ Please do not remove and change any accessories of the machine.
- ◆◆ It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆ It forbids others to change and maintain the machine except the professional people.
- ◆◆ Please contact us if there is any trouble during operation.

**This product does not belong to medical,
please don't use for disease treatment!**

Parts Introduction



- 1.Handle base
- 2.Touch screen
- 3.Service brake
- 4.Beating work handpiece jack
- 5.Service brake jack
- 6.General power switch
- 7.Power supply line switch

Accessories



Beating work handpiece



1# Applicator head



2# Applicator head



3# Applicator head



4# Applicator head



5# Applicator head

The 5# Applicator head needs to be installed with a blue silicone cap, after installation as shown.



Other types of applicator head do not need silicone cap.

Any part of instrument made of PMMA material (for example: accessories shelf, panel, transparent handles, etc.). Do not use alcohol or corrosive solvent to wipe clean.

Basal Operation

1.Put accessories well and connect them to the host.
(1)Insert the plug of Beating work handpiece into [4] output jack.

(2)Insert the plug of service brake into [5] output jack.

(3)Choose appropriate applicator head. Hold the handpiece in one hand and unscrew the applicator head from the work handpiece with the other hand (anticlockwise), screw the new applicator head tightly into the work handpiece (clockwise).

2.Insert the plug of power supply line into [7] socket, turn on general power switch [6], put into the system default password “000000”, the system enter interface as (Fig 1) shows.

3.In interface as (Fig 1):

(1)Click on “Basic setting” to select language, set time and date.

NOTICE: “Advanced setting” is only for engineering debugging, not for user.

(2)Click on the corresponding picture where you want to treatment of the screen to enter function interface.

4.This system include function operation parts. This use’s manual use Elbow treatment for example. In interface as (Fig 1), click on “Elbow” to enter interface as (Fig 2) shows.

5.In the interface as (Fig 2):

(1) “Preselection 2500 Times” stands for when the system works to 2500 times every time, it will pause automatically. Please click on any place of the screen according to the hint.


(2)“Current” stands for the working time of the current treatment course . Click on “RESET” to clear the number of current treatment.

(3) “Frequency 1 Hz” stands for the system default work frequency is 1 Hz, click on “<” and “>” to regulate.
The range is : 1Hz ~16Hz.

(4) “Energy 60 mj” stands for the system default work frequency is 60mj, click on “<” and “>” to regulate.



(Fig 2)




when the frequency is 1Hz to 12Hz, the energy rang is : 60mj ~ 185mj. when the frequency is 13Hz or 14Hz, the energy rang is : 60mj ~ 120mj. when the frequency is 15Hz or 16Hz, the energy rang is : 60mj ~ 90mj.

(5)Click on “▶” key, the system is in standby state shooting. Beautician put the service brake under your feet, step on to start working, move away for pause. Click on “||”, the system pause.

6.This machine can click keys on the touchscreen to set the working parameters, it also can use [2]Function selecting knob & Confirming button to set. Turn the knob when you choose the function, the icon “🍀” move to right place, then press knob[2] to confirm.

7.This use’s manual use touchscreen for example.




(Fig 4)



(Fig 5)




(Fig 6)



(Fig 7)



(Fig 8)



(Fig 9)



(Fig 10)




(Fig 11)

Operation

1. Clean the treatment part.
2. According to the treatment part to select appropriate applicator head.
3. Set suitable frequency (the max range is 1Hz~16Hz) and energy (the max range is 60mj~185mj).
4. Operating on. Beautician hold the Beating work handpiece in one hand and place the applicator head on cellulite to stretch marks just ok. The weight of the work handpiece means that it is normally not necessary to apply pressure to the stretch marks treatment.
5. Clean the leftovers with hot towel, use warm towel to clean the work head after treatment.

Cautions:

1. You can use a little gel or other lubricants on the treatment parts, but do not use too more, preventing the lubricants get into the work handpiece.
2. Despite high internal damping as a result of the weight and design of the handpiece, continuous vibrations may cause strain to the user's hand.
3. To avoid heat accumulation in the work head, it is essential to ensure that the air vents at the top and particularly on the base of the work handpiece are not blocked by the hand holding it or anything else.
4. Applicator heads are expendable parts and must be replaced after a certain period functionality.
5. The metal stick inside Beating handpiece must replace when the service life is above two million times.



Contraindications

1. Application to certain tissues: The eyes and the surrounding area, the myocardium, the spinal cord, the gonads, the kidneys and the liver.
2. Blood disorders, coagulation problems or the use of anticoagulants.
3. Polypus in the area of treatment.
4. Thrombosis.
5. Tumor diseases.
6. Polyneuropathy.
7. Growing cartilage in children.
8. Therapy using corticoids.
9. Inapplicable on areas of the body and organs with possible gas content.
10. Inapplicable on areas in proximity to large nerve bundles, blood vessels, the spinal cord and the head

Possible Side Effects

1. Erythema or swelling can temporarily occur in the treated area.
2. Loss of bodily sensation or itching can temporarily occur in the treated area.
3. Skin damage after previous corticoid therapy.
4. Beating application can cause undesirable heart activity.

Touch Screen Calibration

1. Before turn on the system, please click on the screen, don't move the finger, then turn on the Power Switch.
2. The prompts "Err0" is on the screen, "After releasing the key, we enter into the calibration procedure, Operating under the prompts of the screen". Then release your finger now.
3. Operating under the prompts:
 - (1) Press the cross dot in left-top "+".
 - (2) Press the cross dot in right-top "+".
 - (3) Press the cross dot in right-bottom "+".
 - (4) Click on the point (left top).
 - (5) Click on the point (right bottom).
4. Calibration screen is finished, the system enter "input password" interface.

Forbidden Group

1. People with heart disease or high blood pressure, or configured cardiac pacemaker.
2. Patients with acute inflammation, asthma, deep vein thrombosis, thyroncus, cancer.
3. People with hemorrhagic disease, trauma or who is bleeding.
4. Pregnant women.
5. Medical Plastic parts, or parts with artificial in filling inside.
6. People with metal inside the body(include metallic tooth).
7. Patient with skin inflammation or with edema.
8. People with an abnormal immune system.
9. Numb or insensitive to heat.

Technical Parameters


Input Voltage: AC100~240V 50/60Hz □

Output Power: ≤150VA






Fuse: FL5X20\3.15A with buffer(with spare fuse)

Packing List

1.Host	1PC	□
2.Beating work handpiece	1PC	□
3.Service brake	1PC	□
4.Applicator head	1SET(5PCS)	□
5.Silicone cap	10PCS	□
6.Handle base	1PC	□
7.Power supply line	1PC	□
8.User manual	1COPY	□
9.Hexagonal wrench(M3)	1PCS	□
10.Metal stick	1PCS	□



Symbols and Meanings

Warning Symbol	Meaning
 Danger	If used incorrectly will bring the risk of death or serious injury!
 Warning	If used incorrectly will cause possibility ill of death or serious injury!
 Attention	If used incorrectly will cause personnel injury or damage to the goods!
 Danger	Device exists high pressure risk of high voltage output terminals!
 Warning	If used incorrectly will cause personnel injury!

*Things damage refers to housing, property, livestock and pet damage.

Contraindications

1. The following symptoms were forbidden to use

- (1) Patients who used by Heart rate regulator, cardiac pacemaker, ECS;
- (2) The use of artificial bone; Patients with silicosis;
- (3) Drunk or high fever patients;
- (4) Patients with higher sensitivity mental disease and epilepsy patients
- (5) Patients with no independent expression consciousness and can not take care of themselves.
- (6) Patients with metal materials in body (including metal teeth), pregnant women, cachexia, malignant tumor.
- (7) Severe hypertension and severe patients with cardiovascular and cerebrovascular disease.
- (8) Patients with renal failure.
- (9) Patients with severe infection



2. The following please use under the guidance of professionals

- (1) Women during menstruation;
- (2) Children;
- (3) Hemorrhagic disease or trauma, inflammation, skin disease, skin infection;
- (4) Patients with feeling numb, no response to cold and heat;
- (5) The immune system has defects, scar formation history, physical disease;
- (6) The people who has taken face-lifting or her body has artificial thing.



Maintenance

1. Appliances and medical equipment close to the machine may be resulting in operational failure due to interference.
2. Prohibit use of the equipment in the following environment
 - Hot and humid environment
 - Dusty environment
 - Environmental with much smoke or steam
3. If there comes a failure, please shut down the power immediately, and send it to the distributor or manufacturer for maintenance. Do not repair by yourself.

4. Exclusion of liability.

The Company is not responsible for the following circumstances:

- (1) If fire, earthquake, war occurred;
- (2) The fault of the user, operator error;
- (3) Damage caused by using other method (not included in the use manual).
- (4) Operate without the manual's guide.
5. Do not make the machine close to the wall, keeping a 30CM space around for cooling.

Safety and Environmental Considerations

Toxic and hazardous substances reference list Electronic information product pollution control and management measures Table of contents of toxic and hazardous substances or elements

Part Name	Toxic and Hazardous Substance or Element						Remark
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr(VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)	
Plastic Board	○	○	○	○	○	○	
Printed Circuit Board	×	○	○	○	○	○	
Power Supply	×	○	○	○	○	○	
Accessory	×	○	○	○	○	○	

Declaration:

- The content of the restricted substances in all homogeneous materials of this part is below the limits set by the SJ/T11363-2006 standard.
 - ×
- The content of a homogeneous material in the part of the material is beyond the SJ/T11363-2006 standard requirements. The user's normal use of this product within the time limit of environmental protection, these substances or elements will not occur leakage or sudden change, will not cause the personal, property damage. For such substances or elements, the user shall not dispose of it yourself, please according to government regulations to the relevant departments of the government designated recycling.