

# BOOM DORMED PLUS



## Gynecological examination chair incorporates four electrical movements as follow:

Electric movement of High-low from 60cm to 90cm  
Electric movement of Back section from 0° to 65°  
Electric movement of Seat section from 0° to 25°  
Electric movement of Trend – Antitrend at preselected angle from -12° to +12°

## Bedding Platform:

Divided in two sections made with integrated upholstery mattress washable and disinfected, in a choice of colors

## Base:

Metal heavy construction painted with electrostatic paint.

## Dimensions:

Back section: 86x56cm  
Seat section: 35x56cm  
Leg section: 45x56cm  
Height adjustment: 60-90cm  
Safety load: 185 kg

## Includes the following accessories:

- Pair of swiveling and height adjustable knee rest (Goepel type) with the corresponding clamps.
- Head cushion removable
- Inox basin detachable
- Inox basin holder detachable for disinfection purposes with locking device
- Pair of metallic handles for the patient with plastic edge
- Paper roll holder
- Plastic cover at the back, seat and back section
- Hand control for four movements

