



High Quality Linear Wireless Probe

- Standard Configuration:
- Main host(Prote) 1 Unit
- Internal battery and charger cable 1 Set
- Options:
- Tablet, Phone
- Wireless charge pad

Specifications:

- High quality:192 elements /64channels
- Scanning system: linear array scan:
- Frequency:10/14MHz(PL linear array)
- Depth: 100-200mm, Adjustable
- Screen: Smart phone or tablet screen
- Supporting system: Apple MacOSx /Android version
- Display mode: B, B/M/CFD/D
- Frame rate: 12f/s
- Image gray scaler:256 level
- Image Storage: Using intelligent terminal customers to use (mobile phones, Tablet PC) storage devices
- Measure: distance, area, obstetrics and other
- Power: by built-in battery
- Power consumption: 10W(unfreeze)/4W(freeze)
- Battery working time: 3 hours
- Size: 156mm*60mm*24mmWeight:308gram

